



User Information Sheet

Ronco Item No. **365, 375, 385, 395, 395XX**

Brand Name : **RONCO NITECH®**



Medical
Device
License



Description : Powder Free Nitech® Disposable Glove
 Materials : Unique blend of synthetic materials
 Donning: Powder Free
 Characteristics: Disposable, Non-Sterile, Ambidextrous
 Colour: Blue
 Texture: Smooth

Specification:

| SIZE | | S | M | L | XL | XXL |
|----------------------|-------------------------|---------|------|------|------|-------|
| ITEM NO: | | 365 | 375 | 385 | 395 | 395XX |
| Length: (mm) | ± 10 mm | 240 | 240 | 240 | 240 | 240 |
| Width: (mm) | ± 5 mm | 85 | 95 | 105 | 115 | 125 |
| Palm Thickness (mm) | ± 0.03mm | 0.11 | 0.11 | 0.11 | 0.11 | 0.11 |
| Physical Properties: | Tensile strength (Mpa): | | | | | |
| | Min. Before aging | 11 MPa | | | | |
| | Min. After aging | 11 MPa | | | | |
| | Elongation (%) | | | | | |
| | Min. Before aging | 450% | | | | |
| | Min. After aging | 400% | | | | |
| Pinholes | | AQL 1.5 | | | | |

Packaging :

Size S-XL : 100 Gloves per box, 10 boxes per case. Total 1,000 Gloves per case
 Size XXL : 90 Gloves per box, 10 boxes per case. Total 900 Gloves per case

Product Features :

- Made from a premium blend of synthetic materials unique to RONCO
- Moderately resistant to hydrocarbon oils, fats, detergents and solvents
- Excellent strength and superior elasticity
- Ease of finger movement and reduces hand fatigue
- Beaded cuffs
- Latex Free

Applications :

- | | | | |
|-------------------|------------------|-----------------------|---------------|
| • Food processing | • Medical | • Environmental | • Automotive |
| • Food service | • Laboratory | • Industrial | • Electronics |
| • Meat packing | • Pharmaceutical | • Beauty / Aesthetics | • Printing |
| • Hospitality | • Dental | • Veterinary | |

Storage & Handling :

Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37 °C)

Special Notes :

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).