

User Information Sheet

Item No.: **765S, 765M, 765L, 765XL**

Brand Name :

Blurite 6EC

Description : Disposable Blue Nitrile Examination Glove
 Material : Nitrile
 Donning : Powder Free
 Characteristics : Disposable, Non-Sterile, Ambidextrous
 Colour : Violet
 Texture : Finger Tip Textured



Specification :

Size :	S	M	L	XL
Item Number (Powder Free) :	765S	765M	765L	765XL
Length : (mm) min	290	290	290	290
Width : (mm) ± 10 mm	85	95	105	115
Thickness :	6 mil			
Physical Properties :				
Tensile strength : minimum	14 MPa			
Elongation : minimum	500%			
Pinhole AQL :	AQL 1.5			

Packaging : 50 gloves x 10 dispensers = 500 gloves per case

Product Features :

- Tested against Fentanyl
- Resistant to a wider range of chemicals than latex
- Chemo Tested
- Textured for better grip
- Reduces hand fatigue
- Beaded cuffs
- Latex Free

Applications :

- Emergency Services (EMS)
- Food service
- Hospitality
- Meat packing
- Customs and Border Protection
- Dental care
- Laboratory
- Pharmaceutical
- Veterinary
- Automotive
- Beauty / Aesthetics
- Electronics
- Industrial
- Food processing
- Printing
- Inspection
- Homeland Security
- Correctional Facilities

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

Special Notes :

Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).