



User Information Sheet

Ronco Item No : **762S, 762M, 762L, 762XL, 762XX**

Brand Name : **RONCO OCTOPUS GRIP**

Description : Nitrile Examination Gloves
 Material : Nitrile
 Donning : Powder Free
 Characteristics : Disposable, Non-Sterile, Ambidextrous
 Colour : Black
 Texture : Fully Diamond Embossed



Medical Device License



Specification	Size :	S	M	L	XL	XXL
	Item Number :	762S	762M	762L	762XL	762XX
	Length : (mm) ± 10 mm	240	240	240	240	240
	Width : (mm) ± 10 mm	95	105	115	125	135
	Thickness :	6 mil				
	Physical Properties :					
	Tensile strength : minimum	15 MPa				
	Elongation : minimum	500%				
	AQL (Pinholes)	AQL 2.5				

Packaging : 50 gloves x 10 dispensers = 500 gloves per case

Product Features :

- Raised Diamond Texture
- Excellent strength & grip
- Ease of finger movement
- Superior elasticity
- Reduces hand fatigue
- Resists a wide variety of chemicals
- Latex Free
- Beaded cuffs

Applications :

- Food processing
- Food service
- Hospitality
- Meat packing
- Pharmaceutical
- Veterinary
- Printing
- Cannabis
- Automotive
- Beauty / Aesthetics
- Tatoo
- Electronics
- Industrial
- Laboratory
- Janitorial
- General Purpose

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).