



# User Information Sheet

Ronco Item No : **768S, 768M, 768L, 768XL, 768XX**

Brand Name : **RONCO OCTOPUS GRIP**

Description : Nitrile Examination Gloves  
 Material : Nitrile  
 Donning : Powder Free  
 Characteristics : Disposable, Non-Sterile, Ambidextrous  
 Colour : Orange  
 Texture : Fully Diamond Embossed



Medical Device License



Specification	Size :	S	M	L	XL	XXL
	Item Number :	<b>768S</b>	<b>768M</b>	<b>768L</b>	<b>768XL</b>	<b>768XX</b>
	Length : (mm)	$\pm 10$ mm	240	240	240	240
	Width : (mm)	$\pm 10$ mm	95	105	115	135
	Thickness :	6 mil				
	Physical Properties :					
	Tensile strength : minimum	15 MPa				
	Elongation : minimum	500%				
AQL ( Pinholes )	AQL 2.5					

Packaging : 50 gloves x 10 dispensers = 500 gloves per case

Product Features :

- Raised Diamond Texture
- Reduces hand fatigue
- Excellent strength & grip
- Resists a wide variety of chemicals
- Ease of finger movement
- Latex Free
- Superior elasticity
- Beaded cuffs

Applications :

- Food processing
- Pharmaceutical
- Automotive
- Industrial
- Food service
- Veterinary
- Beauty / Aesthetics
- Laboratory
- Hospitality
- Printing
- Environmental
- Janitorial
- Meat packing
- Cannabis
- Electronics
- General Purpose

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F ( 37°C )

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).