

Ronco Item No : **735X**

Brand Name : **RONCO NE1 (PF)**

Description : 3 mil Nitrile Examination Gloves
 Material : Nitrile
 Donning : Powder Free
 Characteristics : Disposable, Non-Sterile, Ambidextrous
 Colour : Sky Blue
 Texture : Smooth -- Finger Tip Textured



Specification :	Size :	XS	S	M	L	XL
	Item Number :		735S	735M	735L	735XL
	Length : (mm) Minimum		230	230	230	230
	Width : (mm) ± 4 mm		85	95	106	115
	Thickness :	3 mil				
	Physical Properties :					
	Tensile strength : minimum	18 MPa				
	Elongation : minimum	500%				
	Pinholes :	AQL 1.5				

Packaging : 100 gloves x 10 dispensers = 1,000 gloves per case

- Product Features :
- Extremely flexible for bare hand feel
 - Light weight yet incredibly strong
 - Sleek, compact packaging helps cut space and cost
 - Resists a variety of chemicals
 - Textured finger tips for good grip

- Applications :
- Medical
 - Dental
 - Veterinary
 - Pharmaceutical
 - Beauty / Aesthetics
 - Laboratory
 - Emergency services
 - Electronics
 - Emergency services
 - Hospitality
 - Food service

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).