

Ronco Item No : **4935XS, 4935S, 4935M, 4935L, 4935XL, 4935XX**

Brand Name :

RONCO CARE

Description : Nitrile Disposable Gloves
 Material : Nitrile
 Donning : Powder Free
 Characteristics : Disposable, Non-Sterile, Ambidextrous
 Colour : Violet
 Texture : Textured Finger Tips



Specification :

Size :		XS	S	M	L	XL	XXL
Item Number :		4935XS	4935S	4935M	4935L	4935XL	4935XX
Length : (mm)	Minimum	230	230	230	230	230	230
Width : (mm)	± 10 mm	70	80	95	110	120	130
Weight: (g)	± 0.5 g	3.2	3.5	3.8	4.1	4.4	4.7
Thickness :		3 mil					

Physical Properties :

Tensile strength (before aging) : minimum	16 MPa
Tensile strength (after aging) : minimum	14 MPa
Elongation (before aging) : minimum	500%
Elongation (after aging) : minimum	400%
Maximum powder level:	2 mg/dm ²

Quality Level :

Major Defects :	AQL 2.5
Minor Defects :	AQL 4.0
AQL (Pinholes):	AQL 4.0

Packaging : 200 gloves x 10 dispensers = 2,000 gloves per case
 180 gloves x 10 dispensers = 1,800 gloves per case (XXL only)

Product Features :

- 100% more gloves per dispenser for environmentally friendly packaging
- Minimized Storage and Transport cost
- Reduced handling costs
- Resists a wide variety of chemicals
- Excellent strength
- Latex free
- Superior elasticity

Applications :

- Food processing
- Food service
- Meat packing
- Veterinary
- Beauty / Aesthetics
- Agricultural
- Janitorial and Sanitation
- Laboratory
- Electronics

Storage & Handling : Store in original packaging in a cool, dry and well-ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)

Special Notes : Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their own determination as to the suitability of the product for their own intended purpose(s).