



focus on quality

# BLURITE™ XPL 958

Nitrile Examination Gloves, 8 mil, 12”  
Extended Cuff, Blue, Powder Free  
10.5 gram (Size: M)



## Available Sizes

Item #	Sizes
958S	Small
958M	Medium
958L	Large
958XL	X-Large
958XX	XX-Large

Packaging:  
50 Gloves x 10 Dispensers = 500 Gloves Per Case  
40 Gloves x 10 Dispensers = 400 Gloves Per Case (For XX-Large Size)

## Features

- RONCO Blurite™ XPL provides an excellent impermeable barrier of protection against bio-hazardous and infectious materials
- Extra long (12 Inches Long) for forearm protection
- Chemical / Solvent Resistance - High
- Flexible ease of movement
- Comfortable for all-day wear
- Disposable, non-sterile, ambidextrous, beaded cuff, fully textured

## Applications

- Medical / Examination
- Automotive
- Electronics
- Chemical Handling
- EMS / Paramedics
- Pharmaceutical
- Laboratory / Clinical
- Fishery / Aquaculture
- Food & Beverage Processing
- Printing

## Storage & Handling



**CAUTION:** Always keep gloves on while working; while removing gloves care must be taken to not touch the outside contaminated part of gloves with bare hands.

In order to maintain proper hygiene standards, never use worn / discarded / contaminated disposable gloves.

## Compliances



Head Office  
70 Planchet Road,  
Concord, Ontario,  
L4K 2C7

☎ 1-905-660-6700  
✉ [contactus@roncosafety.com](mailto:contactus@roncosafety.com)  
🌐 [www.roncosafety.com](http://www.roncosafety.com)



Follow us @RoncoSafety



Although reasonable care has been taken in the preparation of this document, no warranties are extended and no liability is assumed. The information is solely provided as a general guideline for product use and care. It is up to end-users to make their determination as to the suitability of the product for their intended purposes. This document and the information contained herein is the property of Ronco and may not be used or reproduced without permission.